



SSBP4 rabbit pAb

Catalog No	YP-Ab-09401
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SSBP4
Immunogen	Synthesized peptide derived from human SSBP4 AA range: 213-263
Specificity	This antibody detects endogenous levels of human SSBP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SSBP4;Single-stranded DNA-binding protein 4
Observed Band	39-42 kDa
Calculated Molecular Weight	39 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

