



BARX2 rabbit pAb

Catalog No	YP-Ab-09402
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BARX2
Immunogen	Synthesized peptide derived from human BARX2 AA range: 167-217
Specificity	This antibody detects endogenous levels of human BARX2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BARX2;Homeobox protein BarH-like 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	Highly expressed in adult salivary gland and at much lower levels in mammary gland, kidney and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
