



CLC9A rabbit pAb

Catalog No	YP-Ab-09403
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLEC9A UNQ9341/PRO34046
Immunogen	Synthesized peptide derived from human CLC9A AA range: 33-83
Specificity	This antibody detects endogenous levels of human CLC9A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C type lectin domain family 9;member A;C-type lectin domain family 9 member A;CLC9A_HUMAN;CLEC9A;DNDR1;HEEE9341;PRO3446;UNQ9341;CLC9A;CLEC9A UNQ9341/PRO34046
Observed Band	27kD
Calculated Molecular Weight	25~45kD; 27kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	In peripheral blood highly restricted on the surface of BDCA31(+) dendritic cells and on a small subset of CD14(+) and CD16(-) monocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
