



F198A rabbit pAb

Catalog No	YP-Ab-09404
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM198A C3orf41
Immunogen	Synthesized peptide derived from human F198A AA range: 486-536
Specificity	This antibody detects endogenous levels of human F198A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F198A;FAM198A C3orf41;Protein FAM198A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Endoplasmic reticulum . Golgi apparatus . Membrane, caveola . Requires caveolae biogenesis to be secreted from the endoplasmic reticulum going through the Golgi apparatus where is post-translationally processed to the mature form. .
Tissue Specificity	Expressed in skin, lung and colon (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images