



LHPL3 rabbit pAb

Catalog No	YP-Ab-09410
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LHFPL3 LHFPL4
Immunogen	Synthesized peptide derived from human LHPL3 AA range: 154-204
Specificity	This antibody detects endogenous levels of human LHPL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LHPL3;LHFPL3 LHFPL4;Lipoma HMGIC fusion partner-like 3 protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images