



MEX3A rabbit pAb

Catalog No	YP-Ab-09415
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MEX3A RKHD4
Immunogen	Synthesized peptide derived from human MEX3A AA range: 48-98
Specificity	This antibody detects endogenous levels of human MEX3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RKHD4;MEX 3A;MEX3A;MEX3A RKHD4;RNA-binding protein MEX3A (RING finger and KH domain-containing protein 4)
Calculated Molecular Weight	54 kDa
Cell Pathway	Cytoplasm . Nucleus . Cytoplasm, P-body . Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway.
Tissue Specificity	Highest levels found in fetal brain and testis. Detected also in thymus, salivary gland and uterus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

