



CL054 rabbit pAb

Catalog No	YP-Ab-09416
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C12orf54 HSD-29 HSD-30 HSD29 HSD30
Immunogen	Synthesized peptide derived from human CL054 AA range: 60-110
Specificity	This antibody detects endogenous levels of human CL054
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CL054;C12orf54 HSD-29 HSD-30 HSD29 HSD30;Uncharacterized protein C12orf54
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
