



FRG2B rabbit pAb

Catalog No	YP-Ab-09417
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FRG2B
Immunogen	Synthesized peptide derived from human FRG2B AA range: 196-246
Specificity	This antibody detects endogenous levels of human FRG2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FRG2B;Protein FRG2-like-1 (FSHD region gene 2 protein family member B) (HSA10-FRG2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
