



LTBP2 rabbit pAb

Catalog No	YP-Ab-09418
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LTBP2 C14orf141 LTBP3
Immunogen	Synthesized peptide derived from human LTBP2 AA range: 1080-1130
Specificity	This antibody detects endogenous levels of human LTBP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf141;LTBP3;Latent-transforming growth factor beta-binding protein 2; LTBP 2;LTBP-2;LTBP2;LTBP2 C14orf141 LTBP3;Latent-transforming growth factor beta-binding protein 2 (LTBP-2)
Observed Band	260 kDa
Calculated Molecular Weight	1821 aa, 195 kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed in the aorta (at protein level). Expressed in lung, weakly expressed in heart, placenta, liver and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
