



RNF31 rabbit pAb

Catalog No	YP-Ab-09419
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF31 ZIBRA
Immunogen	Synthesized peptide derived from human RNF31 AA range: 833-883
Specificity	This antibody detects endogenous levels of human RNF31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ10111;FLJ23501;HOIL 1 interacting protein;HOIP;ring finger protein 31; RNF31;ZIBRA;RNF31 ZIBRA;E3 ubiquitin-protein ligase RNF31 (EC 6.3.2.-) (HOIL-1-interacting protein) (HOIP) (RING finger protein 31) (Zinc in-between-RING-finger ubiquitin-associated domain protein)
Observed Band	100 kDa
Calculated Molecular Weight	1072 aa, 120 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in both normal and transformed breast epithelial cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
