



ZBED1 rabbit pAb

Catalog No	YP-Ab-09424
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZBED1 ALTE DREF KIAA0785 TRAMP
Immunogen	Synthesized peptide derived from human ZBED1 AA range: 329-379
Specificity	This antibody detects endogenous levels of human ZBED1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALTE;DREF;dREF homolog;hDREF;KIAA785;TRAMP;ZBED1;ZBED1 ALTE DREF KIAA0785 TRAMP;Zinc finger BED domain-containing protein 1 (Putative Ac-like transposable element) (dREF homolog)
Observed Band	
Calculated Molecular Weight	78 kDa
Cell Pathway	Nucleus, PML body . Nucleus . In granular structures. .; Nucleus, PML body. (Microbial infection) Upon interaction with human adenovirus early E1A protein, the protein is redistributed to the peripheral areas of PML bodies. .
Tissue Specificity	Ubiquitously expressed at low levels (PubMed:9887332). Expression is highest in skeletal muscle, heart, spleen and placenta (PubMed:9887332).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images