



CAB45 rabbit pAb

Catalog No	YP-Ab-09426
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SDF4 CAB45 PSEC0034
Immunogen	Synthesized peptide derived from human CAB45 AA range: 63-113
Specificity	This antibody detects endogenous levels of human CAB45
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	45 kDa calcium binding protein;Cab45;SDF 4;SDF4;stromal cell derived factor 4; SDF4 CAB45 PSEC0034;45 kDa calcium-binding protein (Cab45) (Stromal cell-derived factor 4) (SDF-4)
Observed Band	
Calculated Molecular Weight	42 kDa
Cell Pathway	[Isoform 1]: Golgi apparatus lumen .; [Isoform 5]: Cytoplasm . Cell membrane . Cell projection, bleb. Isoform 5 colocalizes with STX3 and STXBP1 isoform 2 at the plasma membrane and cell surface blebs. .
Tissue Specificity	Ubiquitous. Isoform 5 is expressed in pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
