



NECD rabbit pAb

Catalog No	YP-Ab-09431
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NDN
Immunogen	Synthesized peptide derived from human NECD AA range: 101-151
Specificity	This antibody detects endogenous levels of human NECD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HsT16328;NDN;NECD_HUMAN;Necdin;Necdin homolog;Necdin MAGE family member;necdin related protein;PWCR;NECD
Observed Band	
Calculated Molecular Weight	36kDa; 36kD
Cell Pathway	Perikaryon. Nucleus. Neural perikarya, translocates to the nucleus of postmitotic neurons and interacts with the nuclear matrix.
Tissue Specificity	Almost ubiquitous. Detected in fetal brain, lung, liver and kidney; in adult heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, spleen, thymus, prostate, testis, ovary, small intestine and colon. Not detected in peripheral blood leukocytes. In brain, restricted to post-mitotic neurons.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images