



RAG1 rabbit pAb

Catalog No	YP-Ab-09432
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RAG1 RNF74
Immunogen	Synthesized peptide derived from human RAG1 AA range: 894-944
Specificity	This antibody detects endogenous levels of human RAG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RAG 1;RAG1;RNF74;RAG1 RNF74;V(D)J recombination-activating protein 1 (RAG-1) (RING finger protein 74) [Includes: Endonuclease RAG1 (EC 3.1.-.-);E3 ubiquitin-protein ligase RAG1 (EC 6.3.2.-)]
Observed Band	110-120 kDa
Calculated Molecular Weight	119 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Maturing lymphoid cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
