



TECTB rabbit pAb

Catalog No	YP-Ab-09433
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TECTB
Immunogen	Synthesized peptide derived from human TECTB AA range: 200-250
Specificity	This antibody detects endogenous levels of human TECTB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TECTB;Beta-tectorin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor ; Extracellular side . Secreted, extracellular space, extracellular matrix. Found in the non-collagenous matrix of the tectorial membrane. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
