



KRA41 rabbit pAb

Catalog No	YP-Ab-09436
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KRTAP4-1 KAP4.10 KRTAP4-10 KRTAP4.1 KRTAP4.10
Immunogen	Synthesized peptide derived from human KRA41 AA range: 68-118
Specificity	This antibody detects endogenous levels of human KRA41
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KRA41;KRTAP4-1 KAP4.10 KRTAP4-10 KRTAP4.1 KRTAP4.10;Keratin-associated protein 4-1 (Keratin-associated protein 4-10) (Keratin-associated protein 4.1) (Keratin-associated protein 4.10) (Ultrahigh sulfur keratin-associated protein 4.10)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in the hair follicles.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

