



MKRN3 rabbit pAb

Catalog No	YP-Ab-09440
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MKRN3 D15S9 RNF63 ZNF127
Immunogen	Synthesized peptide derived from human MKRN3 AA range: 237-287
Specificity	This antibody detects endogenous levels of human MKRN3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D15S9;makorin ring finger protein 3;MKRN3;RING finger protein 63;RNF63;ZFP127;Zinc finger protein 127;ZNF127;MKRN3 D15S9 RNF63 ZNF127; Probable E3 ubiquitin-protein ligase makorin-3 (EC 6.3.2.-) (RING finger protein 63) (Zinc finger protein 127)
Observed Band	56 kDa
Calculated Molecular Weight	507 aa, 56 kDa
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
