



YN4602

FBX30 rabbit pAb

Catalog No	YP-Ab-09442
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO30 FBX30
Immunogen	Synthesized peptide derived from human FBX30 AA range: 212-262
Specificity	This antibody detects endogenous levels of human FBX30
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDC4;F-box and WD-40 domain-containing protein 7;F-box protein FBX30;F-box/WD repeat-containing protein 7;FBW7;FBX30;FBXO30 FBX30;F-box only protein 30
Observed Band	100-110 kDa, 66-75 kDa
Calculated Molecular Weight	66 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images