



BOD1 rabbit pAb

Catalog No	YP-Ab-09443
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BOD1 FAM44B
Immunogen	Synthesized peptide derived from human BOD1 AA range: 112-162
Specificity	This antibody detects endogenous levels of human BOD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BOD1;BOD1 FAM44B;Biorientation of chromosomes in cell division protein 1 (Biorientation defective protein 1) (Protein FAM44B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Chromosome, centromere, kinetochore . Localizes at the centrosomes throughout the cell cycle, only dissociating during cytokinesis. Localizes at the kinetochore from prometaphase until anaphase.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images