



AT5EL rabbit pAb

Catalog No	YP-Ab-09445
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ATP5EP2
Immunogen	Synthesized peptide derived from human AT5EL AA range: 1-51
Specificity	This antibody detects endogenous levels of human AT5EL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AT5EL;ATP5EP2;ATP synthase subunit epsilon-like protein;mitochondrial
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
