



YN4606

bQCR10 rabbit pA

Catalog No	YP-Ab-09446
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UQCR11 UQCR
Immunogen	Synthesized peptide derived from human QCR10 AA range: 4-54
Specificity	This antibody detects endogenous levels of human QCR10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UQCR11;Complex III subunit 10;Complex III subunit XI;Cytochrome b c1 complex subunit 10;Cytochrome b-c1 complex subunit 10;QCR10;UQCR11 UQCR;Cytochrome b-c1 complex subunit 10 (Complex III subunit 10) (Complex III subunit XI) (Ubiquinol-cytochrome c reductase complex 6.4 kDa protein)
Observed Band	6 kDa
Calculated Molecular Weight	6 kDa
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images