



DNJC8 rabbit pAb

Catalog No	YP-Ab-09447
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC8 SPF31 HSPC315 HSPC331
Immunogen	Synthesized peptide derived from human DNJC8 AA range: 128-178
Specificity	This antibody detects endogenous levels of human DNJC8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNJC8;DNAJC8 SPF31 HSPC315 HSPC331;DnaJ homolog subfamily C member 8 (Splicing protein spf31)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
