



ANS4B rabbit pAb

Catalog No	YP-Ab-09449
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ANKS4B HARP
Immunogen	Synthesized peptide derived from human ANS4B AA range: 358-408
Specificity	This antibody detects endogenous levels of human ANS4B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANS4B;ANKS4B HARP;Ankyrin repeat and SAM domain-containing protein 4B (Harmonin-interacting ankyrin repeat-containing protein) (Harp)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, microvillus . Localizes at the tip of microvilli (PubMed:26812018). May associate with endoplasmic reticulum membranes (By similarity). .
Tissue Specificity	Expressed in kidney and small intestine.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
