



MEST rabbit pAb

Catalog No	YP-Ab-09452
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MEST PEG1
Immunogen	Synthesized peptide derived from human MEST AA range: 187-237
Specificity	This antibody detects endogenous levels of human MEST
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MEST;PEG1;MEST PEG1;Mesoderm-specific transcript homolog protein (EC 3.-.-.) (Paternally-expressed gene 1 protein)
Observed Band	34 kDa
Calculated Molecular Weight	38 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Highly expressed in hydatidiform moles, but barely expressed in dermoid cysts. Biallelic expression is detected in blood lymphocytes. Seems to imprinted in an isoform-specific manner rather than in a tissue-specific manner in lymphocytes. Isoform 1 is expressed only from the paternal allele. Isoform 2 is expressed from both the paternal allele and the maternal allele.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images