



FA156 rabbit pAb

Catalog No	YP-Ab-09454
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM156A TMEM29 PP12994 PRO0659; FAM156B TMEM29B
Immunogen	Synthesized peptide derived from human FA156 AA range: 6-56
Specificity	This antibody detects endogenous levels of human FA156
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA156;FAM156A TMEM29 PP12994 PRO0659;FAM156B TMEM29B;Protein FAM156A/FAM156B (Transmembrane protein 29)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
