



MBD4 rabbit pAb

Catalog No	YP-Ab-09457
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MBD4 MED1
Immunogen	Synthesized peptide derived from human MBD4 AA range: 421-471
Specificity	This antibody detects endogenous levels of human MBD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.2.2.-;MED1;Methyl-CpG-binding domain protein 4;Methyl-CpG-binding endonuclease 1;Methyl-CpG-binding protein MBD4;MBD4;MBD4 MED1;Methyl-CpG-binding domain protein 4 (EC 3.2.2.-) (Methyl-CpG-binding endonuclease 1) (Methyl-CpG-binding protein MBD4) (Mismatch-specific DNA N-glycosylase)
Observed Band	60-70 kDa
Calculated Molecular Weight	66 kDa
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
