



NUD16 rabbit pAb

Catalog No	YP-Ab-09460
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NUDT16
Immunogen	Synthesized peptide derived from human NUD16 AA range: 2-52
Specificity	This antibody detects endogenous levels of human NUD16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ31265;FLJ3434;FLJ36248;Nucleoside diphosphate-linked moiety X motif 16; NUD16_HUMAN;Nudix (nucleoside diphosphate linked moiety X) type motif 16; Nudix motif 16;NUDT 16;NUDT16;U8 snoRNA-binding protein H29K;U8 snoRNA-decapping enzyme;NUD16;U8 snoRNA-decapping enzyme (EC 3.6.1.62) (Nucleoside diphosphate-linked moiety X motif 16) (Nudix motif 16) (U8 snoRNA-binding protein H29K) (m7GpppN-mRNA hydrolase)
Observed Band	
Calculated Molecular Weight	21 kDa; 21kD
Cell Pathway	Nucleus . Nucleus, nucleoplasm . Nucleus, nucleolus . Cytoplasm . Localized predominantly in the cytoplasm (PubMed:21070968). Localized in nucleolus, and in a minor proportion in distinct foci in the nucleoplasm (By similarity). .
Tissue Specificity	Expressed strongly in lung, kidney, adrenal gland, testis, heart and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images