



TF3C5 rabbit pAb

Catalog No	YP-Ab-09462
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GTF3C5 CDABP0017
Immunogen	Synthesized peptide derived from human TF3C5 AA range: 458-508
Specificity	This antibody detects endogenous levels of human TF3C5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TF3C5;GTF3C5 CDABP0017;General transcription factor 3C polypeptide 5 (TF3C-epsilon) (Transcription factor IIIC 63 kDa subunit) (TFIIIC 63 kDa subunit) (TFIIIC63) (Transcription factor IIIC subunit epsilon)
Observed Band	
Calculated Molecular Weight	6kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
