



HRSL2 rabbit pAb

Catalog No	YP-Ab-09463
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HRASLS2
Immunogen	Synthesized peptide derived from human HRSL2 AA range: 1-51
Specificity	This antibody detects endogenous levels of human HRSL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HRSL2;HRASLS2;HRAS-like suppressor 2 (EC 2.3.1.-) (EC 3.1.1.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Membrane ; Single-pass membrane protein . Exhibits a granular pattern in the cytoplasm with preferential perinuclear localization. .
Tissue Specificity	Expressed in liver, kidney, small intestine testis and colon (PubMed:19615464). Undetectable in testis, placenta, salivary gland and fetal brain (PubMed:18163183).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

