



NUD22 rabbit pAb

Catalog No	YP-Ab-09465
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NUDT22 PP11246
Immunogen	Synthesized peptide derived from human NUD22 AA range: 44-94
Specificity	This antibody detects endogenous levels of human NUD22
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NUD22;NUDT22 PP11246;Nucleoside diphosphate-linked moiety X motif 22 (Nudix motif 22)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
