



INSM1 rabbit pAb

Catalog No	YP-Ab-09466
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	INSM1 IA1
Immunogen	Synthesized peptide derived from human INSM1 AA range: 61-111
Specificity	This antibody detects endogenous levels of human INSM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IA-1;IA 1;IA1;INSM 1;insulinoma associated 1;INSM1;INSM1 IA1;Insulinoma-associated protein 1 (Zinc finger protein IA-1)
Observed Band	53-60 kDa
Calculated Molecular Weight	53 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in pancreatic duct cells. Expressed in several tumor cell lines of neuroendocrine origin including pheochromocytoma, medullary thyroid carcinoma, insulinoma, medulloblastoma, retinoblastoma, pheochromacytoma, medullary thyroid carcinoma and small cell lung carcinoma.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images