



RCN1 rabbit pAb

Catalog No	YP-Ab-09467
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RCN1 RCN
Immunogen	Synthesized peptide derived from human RCN1 AA range: 136-186
Specificity	This antibody detects endogenous levels of human RCN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADAPT78;Calcipressin 1;CSP1;DSC1;DSCR1;MCIP1;RCAN1;RCN1;regulator of calcineurin 1;RCN1 RCN;Reticulocalbin-1
Observed Band	25-29 kDa
Calculated Molecular Weight	28 kDa
Cell Pathway	Endoplasmic reticulum lumen.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
