



MFA3L rabbit pAb

Catalog No	YP-Ab-09469
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MFAP3L KIAA0626 HSD-39 HSD39
Immunogen	Synthesized peptide derived from human MFA3L AA range: 151-201
Specificity	This antibody detects endogenous levels of human MFA3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MFA3L;MFAP3L KIAA0626 HSD-39 HSD39;Microfibrillar-associated protein 3-like (Testis development protein NYD-SP9)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Nucleus . Cytoplasm . Mainly localized in the nucleus (PubMed:24735981). .
Tissue Specificity	Highly expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
