



PRAM rabbit pAb

Catalog No	YP-Ab-09470
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PRAM1
Immunogen	Synthesized peptide derived from human PRAM AA range: 310-360
Specificity	This antibody detects endogenous levels of human PRAM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PML RAR;PRAM 1;PRAM1;PRAM;PML-RARA-regulated adapter molecule 1 (PRAM) (PRAM-1)
Observed Band	70 kDa
Calculated Molecular Weight	718 aa, 79 kDa
Cell Pathway	
Tissue Specificity	Expressed in peripheral blood leukocytes and bone marrow. Expressed in monocytes, and to a lesser extent in granulocytes and lymphocytes. Not expressed in non hematopoietic tissues except in lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
