



EBPL rabbit pAb

Catalog No	YP-Ab-09471
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EBPL EBRP ERP
Immunogen	Synthesized peptide derived from human EBPL AA range: 10-60
Specificity	This antibody detects endogenous levels of human EBPL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EBPL;EBRP;emopamil binding protein like;ERP;EBPL EBRP ERP;Emopamil-binding protein-like (Emopamil-binding-related protein)
Observed Band	
Calculated Molecular Weight	26 aa, 23 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed with highest levels in liver, lung and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
