



ASB12 rabbit pAb

Catalog No	YP-Ab-09473
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ASB12
Immunogen	Synthesized peptide derived from human ASB12 AA range: 113-163
Specificity	This antibody detects endogenous levels of human ASB12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ankyrin repeat and SOCS box protein 12;ASB-12;ASB12;Ankyrin repeat and SOCS box protein 12 (ASB-12)
Observed Band	34,39 kDa
Calculated Molecular Weight	309 aa, 34 kDa/318 aa, 35 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
