



TPIS rabbit pAb

Catalog No	YP-Ab-09475
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPI1 TPI
Immunogen	Synthesized peptide derived from human TPIS AA range: 76-126
Specificity	This antibody detects endogenous levels of human TPIS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPAG1;sperm associated antigen 1;SP75;Infertility-related sperm protein Spag-1;HSD-3.8;TPIS;TPI1 TPI;Triosephosphate isomerase (TIM) (EC 5.3.1.1) (Triose-phosphate isomerase)
Observed Band	104 kDa
Calculated Molecular Weight	104 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
