



H2B1J rabbit pAb

Catalog No	YP-Ab-09476
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HIST1H2BJ H2BFR
Immunogen	Synthesized peptide derived from human H2B1J AA range: 58-108
Specificity	This antibody detects endogenous levels of human H2B1J
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	dJ221C16.8;H2B histone family;member A;H2B/a;H2B/g;H2B/h;H2B/k;H2B/l;H2B1C_HUMAN;H2BFA;HIST1H2BI;Histone 1;H2bg;Histone cluster 1 H2B family member g;Histone cluster 1;Histone H2B type 1-C/E/F/G/I;Histone H2B.1 A;Histone H2B.a;Histone H2B.g;Histone H2B.h;Histone H2B.k;Histone H2B.l;H2B histone family member R;H2B/r;H2B1J_HUMAN;H2BFR;HIST1H2BJ;Histone 1 H2bj;Histone cluster 1 H2bj;Histone H2B type 1-J;Histone H2B.1;Histone H2B.r;GL15;H2B;H2B histone family member Q;H2B-GL15;H2B.1;H2B/q;H2B2E_HUMAN;H2BFQ;H2BGL15;H2BQ;HIST2H2BE;Histone 2 H2be;histone cluster 2;H2be;histone H2B GL15;histone H2B type 2 E;Histone H2B type 2-E;histone H2B type 2E;Histone H2B-GL15;Histone H2B.q;histone H2BGL15;H2B K;H2B type 12;H2B/b;H2B/c;H2B/d;H2B/e;H2B/f;H2B/j;H2B/n;H2B/s;HIRA-interacting protein 1;HIRA-interacting protein 2;Histone H2B type 1-B;Histone H2B type 1-D;Histone H2B type 1-H;Histone H2B type 1-K;Histone H2B type 1-L;Histone H2B type 1-M;Histone H2B type 1-N;Histone H2B type 1-O;Histone H2B type 2-F;Histone H2B type 3-B;Histone H2B type F-S;Histone H2B.1 B;Histone H2B.2;Histone H2B.b;Histone H2B.c;Histone H2B.d;Histone H2B.e;Histone H2B.f;Histone H2B.j;Histone H2B.n;Histone H2B.s;H2B1J;HIST1H2BJ H2BFR;Histone H2B type 1-J (Histone H2B.1) (Histone H2B.r) (H2B/r)
Observed Band	
Calculated Molecular Weight	14kDa; 14kD



Cell Pathway	Nucleus. Chromosome.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
