



SNG1 rabbit pAb

Catalog No	YP-Ab-09477
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SYNGR1
Immunogen	Synthesized peptide derived from human SNG1 AA range: 152-202
Specificity	This antibody detects endogenous levels of human SNG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SNG1_HUMAN;Synaptogyrin-1;Syngr1;SNG1
Observed Band	
Calculated Molecular Weight	25kD
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane ; Multi-pass membrane protein . Melanosome . Identified by mass spectrometry in melanosome fractions from stage I to stage IV. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

