



FBX3 rabbit pAb

Catalog No	YP-Ab-09478
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO3 FBX3
Immunogen	Synthesized peptide derived from human FBX3 AA range: 245-295
Specificity	This antibody detects endogenous levels of human FBX3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F box only protein 3;F box protein 3;FBA;FBX3;FBXO3;FBXO3 FBX3;F-box only protein 3
Observed Band	48-50 kDa
Calculated Molecular Weight	471 aa, 55 kDa
Cell Pathway	Nucleus. Colocalizes with PML at the peripheries of nuclear bodies.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
