



TBCA rabbit pAb

Catalog No	YP-Ab-09480
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TBCA
Immunogen	Synthesized peptide derived from human TBCA AA range: 45-95
Specificity	This antibody detects endogenous levels of human TBCA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CFA;TBCA;TCP1 chaperonin cofactor A;tubulin folding cofactor A;Tubulin specific chaperone A;Tubulin-specific chaperone A (TCP1-chaperonin cofactor A) (Tubulin-folding cofactor A) (CFA)
Observed Band	13 kDa
Calculated Molecular Weight	108 aa, 13 kDa
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
