



TMOD4 rabbit pAb

Catalog No	YP-Ab-09484
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMOD4
Immunogen	Synthesized peptide derived from human TMOD4 AA range: 90-140
Specificity	This antibody detects endogenous levels of human TMOD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Sk-Tmod;Tropomodulin 4;Tropomodulin-4;TMOD4;Tropomodulin-4 (Skeletal muscle tropomodulin) (Sk-Tmod)
Observed Band	39 kDa
Calculated Molecular Weight	345 aa, 39 kDa
Cell Pathway	Cytoplasm, cytoskeleton . In myofibrils with sarcomeric structure, localizes to the pointed end of actin thin filaments (PubMed:25250574). .
Tissue Specificity	Highly expressed in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
