



AGAL rabbit pAb

Catalog No	YP-Ab-09486
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GLA
Immunogen	Synthesized peptide derived from human AGAL AA range: 138-188
Specificity	This antibody detects endogenous levels of human AGAL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	beta galactosidase 2;BETA GALACTOSIDASE PROTECTIVE PROTEIN;beta-galactosidase 2;beta-galactosidase protective protein;betagalactosidase 2; Carboxypeptidase C;Carboxypeptidase L;carboxypeptidase Y-like kininase; Cathepsin A;Ctsa;deamidase;EC 3.4.16.5;Glactosialidosis;GLB2;Goldberg Syndrome;GSL;lysosomal carboxypeptidase A;Lysosomal protective protein 2 kDa chain;Lysosomal protective protein;Lysosomal protective protein deficiency; NEURAMINIDASE BETA GALACTOSIDASE EXPRESSION;NGBE; Neuraminidase deficiency with beta-galactosidase deficiency;OTTHUMP31778; OTTHUMP31781;PPCA;PPCA deficiency;PPGB;PPGB_HUMAN;Protective protein cathepsin A;Protective protein for beta galactosidase;Protective protein for beta-galactosidase;Protective protein/cathepsin A deficiency;urinary kininase; AGAL;GLA;Alpha-galactosidase A (EC 3.2.1.22) (Alpha-D-galactosidase A) (Alpha-D-galactoside galactohydrolase) (Melibiase) (Agalsidase)
Calculated Molecular Weight	54kDa; 54kD
Cell Pathway	Lysosome.
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

