



PRUNE rabbit pAb

Catalog No	YP-Ab-09488
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRUNE
Immunogen	Synthesized peptide derived from human PRUNE AA range: 48-98
Specificity	This antibody detects endogenous levels of human PRUNE
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRES-17;DRES17;Drosophila-related expressed sequence 17;hPrune;HTcD37; Protein prune homolog;PRUNE;Prune homolog;PRUNE like protein; PRUNE_HUMAN;TcD37 homolog;Protein prune homolog (hPrune) (EC 3.6.1.1) (Drosophila-related expressed sequence 17) (DRES-17) (DRES17) (HTcD37)
Observed Band	
Calculated Molecular Weight	5kD
Cell Pathway	Cytoplasm. Nucleus. Cell junction, focal adhesion. In some transfected cells a nuclear staining is also observed.
Tissue Specificity	Ubiquitously expressed. Seems to be overexpressed in aggressive sarcoma subtypes, such as leiomyosarcomas and malignant fibrous histiocyctomas (MFH) as well as in the less malignant liposarcomas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images