



DYH8 rabbit pAb

Catalog No	YP-Ab-09493
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DNAH8
Immunogen	Synthesized peptide derived from human DYH8 AA range: 2086-2136
Specificity	This antibody detects endogenous levels of human DYH8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DYH8;DNAH8;Dynein heavy chain 8;axonemal (Axonemal beta dynein heavy chain 8) (Ciliary dynein heavy chain 8)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, flagellum axoneme . Cytoplasm . Detected in sperm tail, with almost exclusive localization to the principal piece. Also detected in the cytoplasm of primary spermatocytes. .
Tissue Specificity	Expressed in spermatozoa (at protein level). Not detected in airway epithelial cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images