



PEPD rabbit pAb

Catalog No	YP-Ab-09494
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PEPD PRD
Immunogen	Synthesized peptide derived from human PEPD AA range: 138-188
Specificity	This antibody detects endogenous levels of human PEPD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.13.9;PRD;PROLIDASE;Xaa-Pro dipeptidase;X-Pro dipeptidase;PEPD; PEPD PRD;Xaa-Pro dipeptidase (X-Pro dipeptidase) (EC 3.4.13.9) (Imidodipeptidase) (Peptidase D) (Proline dipeptidase) (Prolidase)
Observed Band	55 kDa
Calculated Molecular Weight	493 aa, 55 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
