



RGCC rabbit pAb

Catalog No	YP-Ab-09495
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RGCC C13orf15 RGC32
Immunogen	Synthesized peptide derived from human RGCC AA range: 5-55
Specificity	This antibody detects endogenous levels of human RGCC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C13orf15;bA157L14.2;RGCC;RGC-32;RGC 32;RGCC C13orf15 RGC32; Regulator of cell cycle RGCC (Response gene to complement 32 protein) (RGC-32)
Calculated Molecular Weight	15 kDa
Cell Pathway	Cytoplasm. Nucleus. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasmic in unstimulated cells. Nuclear after activation by complement. Associated with the centrosome during prometaphase and metaphase.
Tissue Specificity	Detected in brain, heart and liver (at protein level). Highly expressed in liver, skeletal muscle, kidney and pancreas. Detected at lower levels in heart, brain and placenta. Detected in aorta endothelial cells. Overexpressed in colon, breast, prostate, bladder, lung, and ovarian cancer tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images