



RNF17 rabbit pAb

Catalog No	YP-Ab-09496
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF17 TDRD4
Immunogen	Synthesized peptide derived from human RNF17 AA range: 926-976
Specificity	This antibody detects endogenous levels of human RNF17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Mmp 2;ring finger protein 17;RNF17;SPATA23;TDRD4;RNF17 TDRD4;RING finger protein 17 (Tudor domain-containing protein 4)
Observed Band	184 kDa
Calculated Molecular Weight	1623 aa, 185 kDa
Cell Pathway	Cytoplasm . Nucleus . Predominantly found in the cytoplasm. Component of a nuage in male germ cells (an electron-dense spherical cytoplasmic body present in late pachytene and diplotene spermatocytes and in elonging spermatids) (By similarity). .
Tissue Specificity	Testis specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images