



RIM3A rabbit pAb

Catalog No	YP-Ab-09497
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RIMBP3 KIAA1666 RIMBP3A
Immunogen	Synthesized peptide derived from human RIM3A AA range: 822-872
Specificity	This antibody detects endogenous levels of human RIM3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RIM3A;RIMBP3 KIAA1666 RIMBP3A;RIMS-binding protein 3A (RIM-BP3.A) (RIMS-binding protein 3.1) (RIM-BP3.1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . In elongating spermatids, localizes to the manchette. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
