



TPRGL rabbit pAb

Catalog No	YP-Ab-09498
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPRG1L FAM79A MOVER
Immunogen	Synthesized peptide derived from human TPRGL AA range: 147-197
Specificity	This antibody detects endogenous levels of human TPRGL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPRGL;TPRG1L FAM79A MOVER;Tumor protein p63-regulated gene 1-like protein (Mossy fiber terminal-associated vertebrate-specific presynaptic protein) (Protein FAM79A)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane ; Peripheral membrane protein . Cell junction, synapse, presynaptic active zone .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
